

<b>Policy #:</b>	701 (PLH-701-04)	<b>Effective Date:</b>	9/30/2004	<b>Reviewed Date:</b>	2/4/2011
<b>Subject:</b>	TRANSFUSION SERVICE SCHEDULE (FROM RECEIPT OF SAMPLE IN TRANSFUSION SERVICE )				
<b>Approved by:</b> Laboratory Executive Director, Ed Hughes (electronic signature)					
<b>Approved by:</b> Laboratory Medical Director, Mark P. Burton, MD (electronic signature)					

## TRANSFUSION SERVICE TEST SCHEDULE (FROM RECEIPT OF SAMPLE IN TRANSFUSION SERVICE)

TRANSFUSION SERVICE TEST	ROUTINE REPORTING INTERVAL	STAT REPORTING INTERVAL	COLLECTION TUBE	SPECIAL INSTRUCTIONS
ABO TYPE	4H	30 MIN	PINK	CHECK HOSPITAL NUMBER
ABO/RH	4H	30 MIN	PINK	CHECK HOSPITAL NUMBER
ANTIBODY SCREEN	4H	45 MIN	PINK	CHECK HOSPITAL NUMBER. SAME TEST AS INDIRECT COOMBS. EXPECTED RANGE:NEGATIVE
ANTIBODY TITER	24H	4H	PINK	CHECK HOSPITAL NUMBER
COLD AGGLUTINATION	6H	1 H	PINK or RED	PERFORMED: 24/7. EDTA TUBE SHOULD BE PLACED IN WARM SAND UNTIL TRANSPORTED TO THE BLOOD BANK.

<b>TRANSFUSION SERVICE TEST</b>	<b>ROUTINE REPORTING INTERVAL</b>	<b>STAT REPORTING INTERVAL</b>	<b>COLLECTION TUBE</b>	<b>SPECIAL INSTRUCTIONS</b>
CROSSMATCH (IS)	1H (current antibody screen sample) 4 H (new sample)	10 MIN (current antibody screen sample, phone order to Blood Bank) 1 H (new sample)	PINK	CHECK HOSPITAL NUMBER, IF PATIENT HAS DESIGNATED DONOR, PLEASE COMMENT WITH ORDER.
DIRECT COOMBS	4H	40 MIN	LAV OR PINK	TEST SHOULD BE PERFORMED WITHIN 24 HOURS OF COLLECTION TIME. EXPECTED RANGE:NEGATIVE
RH TYPE	4H	30 MIN	PINK	CHECK HOSPITAL NUMBER
RHIG WORKUP-BABY	24H	1H (if all samples are in lab)	CORD BLOOD – PINK	MICRO EDTA HEELSTICK IF CORD BLOOD IS NOT AVAILABLE
RHIG-MOM	24H	1H (if all samples are in lab)	PINK	CHECK WITH BLOOD BANK BEFORE DRAWING SPECIMEN
THERAPEUTIC PHLEBOTOMY	24H	N/A	N/A	INCREASED HGB OR HCT REQUIRED
TRANSFUSION REACTION WORKUP	N/A	90 MIN	PINK	SEE PROCEDURE IN NURSING MANUAL
TYPE & SCREEN	4H	45 MIN	PINK	CHECK HOSPITAL NUMBER. NO BLOOD IS CROSSMATCHED WITH THIS ORDER